




Cell-free DNA Analysis Solution

LiquidSCAN

geninus 

© 2021.GENINUS INC.All rights reserved



What is Liquid Biopsy?

A liquid biopsy is a simple and non-invasive alternative to tissue biopsy which enables researchers to profile genetic information of tumor by using DNA extracted from blood.


Utilizing a liquid biopsy, LiquidSCAN® analyzes circulating cell-free DNAs(cfDNAs), the degraded DNA fragments released from various cells to the blood, to detect tumor - derived DNA among cfDNAs.

Due to its non-invasive nature, it can be used for those who are unable or difficult to obtain tumor tissue for diagnosis, therapy selection, treatment monitoring, and detection of minimal residual disease and recurrence.

What is LiquidSCAN®?

LiquidSCAN® is a hybrid capture-based, next-generation sequencing(NGS) assay designed to enable researchers to interrogate genetic variants in cfDNA across cancer related genes. The assay identifies all classes of alterations including single nucleotide variants, indels, copy number alterations, and fusion at a high sensitivity and specificity, which is essential to accurate identification of ctDNA.

Workflow



Streamlined analysis of cfDNA across cancer-related genes. Libraries are constructed with a proprietary assay using molecular barcodes. Illumina NGS is subsequently performed. In addition to our streamlined variant detection analysis, integration for customized analysis is also available upon request.

Types of variants detected by LiquidSCAN®

| Category | Number of genes | SNV and Indels | Fusions | CNV | MSI | TMB |
|---------------------|-----------------|----------------|---------|-----|-----|-----|
| LiquidSCAN® Compact | 44 | ✓ | ✓ | ✓ | ✓ | - |
| LiquidSCAN® IO+ | 243 | ✓ | ✓ | ✓ | ✓ | ✓ |

Features and benefits

Accurate

- Molecular barcoding
- *In silico* digital error suppression algorithm
- Improved library production efficiency

Flexible

- Custom panels
- Module-based analysis service

Comprehensive

- Streamlined assistance from study design to data analysis
- OncoSTATION™ - from sample requisition to data management


Reliable

- College of American Pathologists (CAP) accreditation
- Good Manufacturing Practice (GMP) facility

Application

1. Monitoring

Unlike tissue biopsy, liquid biopsy allows to interrogate clonally divergent and distant lesions without sampling bias thus suits for longitudinal monitoring of patients during treatment. In line with this, LiquidSCAN® is utilized for treatment monitoring, drug resistance monitoring and minimal residual disease monitoring by analyzing ctDNAs at clinical level.



2. Therapy guidance

With Geninus proprietary analysis algorithm, LiquidSCAN® analysis report delivers actionable insights to clinicians. Variants are classified in tiers established by American College of Medical Genetics (ACMG), the targeted drug annotations according to each mutation, and clinical trial information are provided to help clinician to select the therapy to the patients.

SNV/INDEL

| Variants info | | | cfDNA | | gDNA | | Supporting Info | | |
|---------------|--------|------------|--------|-------|--------|-------|-----------------|-------------------------------------------------------------|-------------|
| Tier | Gene | Alteration | VAF(%) | Dep. | VAF(%) | Dep. | Cancer Type | Drug | Sensitivity |
| 1A | PIK3CA | C420R | 13.38 | 2,302 | 0 | 1,828 | Breast Cancer | Alpelisib + Fulvestrant | Sensitive |
| 1B | PIK3CA | C420R | 13.38 | 2,302 | 0 | 1,828 | Breast Cancer | Buparlisib, Serabelisib, Copanlisib, Taselisib+ Fulvestrant | Sensitive |

* Example only, result may differ by cases.

Strength


Improved library production efficiency

Liquid biopsies have a very limited amount of DNA, minimizing the loss between experiments is very important for data analysis. Geninus improved the library production efficiency by 2.7 times or more compared to the standard method of commercial kit by optimizing the reaction conditions of the adaptor ligation step.



Proprietary molecular barcoding technology

To improve data production efficiency, an independent molecular barcoding structure was developed to suppress errors caused by barcode switching that occurred in the traditional molecular barcoding. We have also developed proprietary algorithm that eliminate cfDNA specific error pattern.



Possessing and customization of panels by cancer type

By using in-house developed database, cancer specific panel such as lymphoma, renal cell carcinoma, prostate cancer, follicular lymphoma and cancer specific noise filtering algorithm is developed.

Based on this panel development capability, we can design a unique panel that meets the needs of customers.



Clinical report : Bridging data and clinicians



- To respond to clinical needs, variants are classified in tiers established by ACMG and contents are designed by careful arrangement for actual usage and flexibility.
- From general information of patient to the drug annotations, LiquidSCAN® report reflects our duty to help clinicians make the most appropriate decision by having the confidence to leverage NGS technology.

OncoSTATION™ : Beyond data, towards clinics




OncoSTATION™ is a software having its competency in supporting Laboratory information management system (LIMS), bioinformatics analysis and visualization UI provided to the customers as a package. Its automated system enables straightforward data analysis using proprietary algorithm. OncoSTATION™ provides a user-friendly environment for pathologists and clinicians by implementing it in sync with Electronic Medical Record (EMR) so that pathologists and clinicians can conveniently apply the test result to the clinical trials and share professional insights with specialists from various fields at the tumor board.

Service specification

Types of variants detected by LiquidSCAN®

| Parameter | | LiquidSCAN Compact | | LiquidSCAN IO+ | |
|----------------------------------|---------------|--------------------------------------------------------------------------------------------------------------------|--------|----------------------------|--------|
| Number of Genes | | 44 | | 243 | |
| Sequencer | | Illumina NextSeq 550Dx or Illumina NextSeq 550 or Illumina NovaSeq 6000 | | | |
| Total run time (Assay to report) | | Fee-for-service at Genius : 4 weeks upon recipient of sample On-premise : 7 days (DNA extraction to sequencing) | | | |
| Kit size | | 96 samples | | | |
| Limit of detection | | 0.5% | | | |
| Detected variant | | SNV / INDEL / FUSION / CNV | | SNV / INDEL / FUSION / CNV | |
| | | MSI | | MSI / TMB | |
| Analytical sensitivity | Variant class | SNV | INDEL | FUSION | CNV |
| | Sensitivity | ≥ 98% | ≥ 98% | ≥ 94% | ≥ 100% |
| | VAF | ≥ 0.5% | ≥ 0.5% | ≥ 0.5% | ≥ 10% |
| Average NGS data output | | 5G | | 20G | |
| Average sequencing depth | | 3000x | | | |

Multiple cancer types and the biomarkers included in LiquidSCAN IO+



Lung

ALK BRAF CDKN2A EGFR KRAS
MAP2K1 NF1 NRAS PIK3CA PTEN
RB1 RET RIT1 SMARCA4 STK11
TP53

Melanoma

ATM BAP1 BRAF CDKN2A MITF
NRAS TERT TP53

Breast

BRCA1 BRCA2 ATM BARD1 BRIP1
CDH1 CHEK2 NF1 PALB2 PTEN
TP53

Colon

PTEN APC CDH1 SMAD4 MLH1
MSH2 MSH6 PMS2 CHEK2 MUTYH
POLD1 POLE STK11

Gastric

CDH1 MLH1 MSH2 MSH6 PMS2
SMAD4

Pancreas

BRCA1 BRCA2 CDKN2A TP53 MLH1
ATM

Sarcoma

TP53 RB1 CDKN2A ATRX CDKN2B
EWSR1 NF1 PTEN MDM2 KMT2D TERT
BRCA2 ARID1B BCOR KDR PIK3CA
CDK4 NTRK1 NOTCH1 MYC ARID1A
NRAS IGF1R DICER1 MET MAP2K4

Liver

TP53 CTNNB1 SMAD2 SMAD4 RB1
CDKN2A CCND1 IDH1 KRAS ARID1A
FGFR2 BAP1 CDKN2B KMT2D ATM
PIK3CA BRAF

Bladder

ARID1A ATM CDKN2A CREBBP
FGFR3 HRAS KDM6A KMT2D
PIK3CA PTEN RAF1 RB1 TP53

Myeloid

KIT FLT3 NPM1 CEBPA RAS WT1
BAALC ERG MN1 DNMT TET2 IDH
ASXL1 PTPN11 CBL

Ovarian

AKT1 BRCA1 BRCA2 CDH1 CTNNB1
MLH1 MSH2 MSH6 PIK3CA PMS2
STK11 TP53

*Representative biomarkers for each cancer type are listed.


Service order process

1) Select service

| | Service | Cancer Type | Number of Genes |
|---|--------------------------------|----------------------|-----------------|
| 1 | LiquidSCAN Compact | Pan-cancer | 44 |
| 2 | LiquidSCAN IO+ | Pan-cancer | 243 |
| 3 | LiquidSCAN Prostate | Prostate | 59 |
| 4 | LiquidSCAN Lymphoma | Lymphoma | 54 |
| 5 | LiquidSCAN RCC | Renal cell carcinoma | 51 |
| 6 | LiquidSCAN Follicular Lymphoma | Lymphoma | 58 |

* Lists of genes for each panel are attached in the Appendix.

2) Sample preparation




Blood

- Two tubes of whole blood (8.5-10ml) are recommended.
(In total, 17~20ml of whole blood are recommended)
- If blood samples are collected in
 - EDTA tubes, plasma samples are needed to be isolated within 6h of collection.
 - Cell-Free DNA™ tubes (Streck inc., Omaha, NE, USA) ships to Geninus in a room temperature storage condition within 2 days (maximum 7 days) after sample collection.



Plasma / buffy coat

- Two tubes of platelet-free plasma (3-5ml per tube) are required
- For matched-normal analysis, one tube of buffy coat samples isolated is required.



DNA extracted from plasma and buffy coat

- Two tubes of 20ng (at least 10ng) cfDNA are desired.
- At least 200ng genomic DNA is required.

3) Shipping

- Cell-Free DNA™ tube : Pack in sealed box and ship at room temperature
- Plasma, buffy coat, DNA : Pack with dry ice in Styrofoam box

📍 5F, KDU Tower, 70, Jeongui-ro, Songpa-gu, Seoul, 05836, Republic of Korea (서울시 송파구 정의로 70(문정동) KDU타워 5층)

☎ 02-6949-6570

Panel and reagent order process

Panel and library preparation kits are available for on-premise analysis.



1) Select panels and reagent

| Product | | | No. of reactions |
|-------------|------------------|-----------------------|------------------|
| LiquidSCAN® | Customized panel | Customized probe | 96 |
| | Library prep kit | EA buffer, enzyme | |
| | | Ligase buffer, enzyme | |
| | | Primer | |
| | | UMI index | |
| | | Blocker solution | |
| | | Binding buffer | |
| | | Wash buffer | |
| | | Purification beads | |

Customized panel and library prep kit can be purchased separately.

2) Analysis software

Analysis software will be setup with OncoSTATION™.

Please contact sales representative for software license.

Contact us

www.kr-genius.com

☎ 02-6949-6570 Fax. 02-6949-6580

✉ info@kr-genius.com

📍 4-5F, KDU Tower, 70, Jeongui-ro, Songpa-gu, Seoul, 05836, Republic of Korea(서울시 송파구 정의로 70(문정동) KDU타워 4-5층)

geninus 